TECHNICAL DATA SHEET

CELOX™ GAUZE

Product Data Sheet
CELOX™ TRAINING GAUZE (5ft)

DESCRIPTION
CELOX™ Training Gauze is a cost saving way to demonstrate wound packing of CELOX™ Rapid Gauze 5ft Z-Fold and CELOX™ Gauze 5ft Z-Fold. It is intended to teach and improve muscle memory for treatment of arterial injuries, road traffic accident and other life-threatening injuries.



PRESENTATION
One 'Z-fold' gauze strip 1.5m long x 7.6cm wide (5ft x 3in).

FEATURES
- Inert moulage trainer for: CELOX™ Rapid Gauze 5ft Z-Fold and CELOX™ Gauze 5ft Z-Fold.
- For simulated training use only.
- Not for medical use.

HOW TO USE
Tightly pack to source of the simulated bleeding, above 'skin level.' Apply a firm compression for three minutes, or until bleeding stops.

INDICATIONS FOR USE
To be used for training purposes only. Not to be used as a medical device on humans to treat major haemorrhage.

ORDERING INFORMATION
FG Code: FG08839030
100 packs per case

STORAGE
CELOX™ products should be stored in dry conditions at ambient temperature.
No special storage conditions are required or indicated on the product labelling.

Shelf life not applicable as this is an unsterile training gauze.

	Sizes
Packet size (approx)	Width – 11.1cm
	Depth – 0.5cm
	Length – 13.2cm
	Weight –32g
Shipper size (approx)	Width – 40cm
	Height – 16cm
	Length – 50cm
	Weight – 3.9kg

CELOX™ MEDICAL
MADE TO SAVE LIVES

For more information visit:
celoxmedical.com

MT-24-200